



# DYNO™

Active lifestyles need active strollers. The Dyno provides positioning for postural support while allowing controlled dynamic movement of the child.

The Dyno folds in one easy step, making life easier for families on the go.

**MAX WEIGHT CAPACITY 80lbs (36kg)**

- 3 WHEEL CONFIGURATION
- LIGHTWEIGHT
- ACTIVATOR™ DYNAMIC SEATING
- CONFIGURABLE SEATING
- 30° HIP ANGLE RANGE
- FOLDS COMPACTLY

The Leggero Dyno is available in 5 bold fabric colors.



Safe



Adjustable



Dynamic



Compact



[www.leggero.us](http://www.leggero.us)

☎ 01869 324600

✉ [simon@stealth-products.co.uk](mailto:simon@stealth-products.co.uk)

# Freedom to go wherever you want... with the positioning options you need.



Complete Seating Options



Dynamic Capabilities



Swiveling Front Wheel



Pneumatic Wheels

**DYNO**

WEIGHT CAPACITY

80 lbs (36 kg)

SEAT WIDTH

14 in (35 cm)

SEAT DEPTH

6-14 in (15-35 cm)

BACK HEIGHT

13-21 in (33-53 cm)

SEAT TO FLOOR

19 in (48 cm)

ADJUSTABLE SEAT TILT

25°-30°

KNEE ACTIVE RANGE  
(OPTIONAL)

10°

SEAT TO BACK ACTIVE RANGE  
(OPTIONAL)

30°

HIP ANGLE ADJUSTMENT

85°, 90°, 105°

LOWER LEG LENGTH

8-14 in (20-35 cm)

CHEST WIDTH RANGE

7-14 in (18-35 cm)

PUSH HANDLE HEIGHT

39-46 in (99-117 cm)

REAR WHEEL DIAMETER

16 in (41 cm)

FRONT WHEEL DIAMETER

10 in (25 cm)

CHAIR WEIGHT

27 lbs (12 kg)

CHAIR WIDTH

26.5 in (67 cm)

FOLDED DIMENSIONS  
(WITH WHEELS)

42 x 17 x 26.5 in (107 x 43 x 67 cm)

FOLDED DIMENSIONS  
(WITHOUT WHEELS)

35 x 17 x 25 in (89 x 43 x 64 cm)



Leggero